
PAβN dihydrochloride

产品编号: D53936

CAS: 100929-99-5

分子式: C₂₅H₃₂Cl₂N₆O₂

纯度: 99.87%

InChi Key: InChiKey=MRUMOHLDMZGMS-IXOXMDGESA-N

Smiles: Cl.Cl.N[C@@H](CC1=CC=CC=C1)C(=O)N[C@@H](CCCNC(N)=N)C(=O)NC1=CC2=CC=CC=C2C=C1

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Antibacterial

溶解性: DMSO:95mg/ml (182.88 mM); H₂O:14.29 mg/mL (27.51 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PAβN dihydrochloride is an inhibitor of efflux pump.