

---

## Fenofibric acid

产品编号: D53938

CAS: 42017-89-0

分子式: C<sub>17</sub>H<sub>15</sub>ClO<sub>4</sub>

纯度: 99.94%

InChi Key: InChIKey=MQOBSOSZFYZQOK-UHFFFAOYSA-N

Smiles: c1c(ccc(c1)C(=O)c1ccc(cc1)OC(C(=O)O)(C)C)Cl

外观: 固体粉末

作用通路: DNA Damage/DNA Repair; Metabolism; Immunology/Inflammation; Neuroscience;  
Proteases/Proteasome

作用靶点: COX; MMP; PPAR

溶解性: DMSO: ≥32 mg/mL (100.39 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Fenofibric Acid is the active form of fenofibrate, a synthetic phenoxy-isobutyric acid derivate with antihyperlipidemic activity.