
Iperoxo

产品编号: D53944

CAS: 247079-84-1

分子式: C₁₀H₁₇N₂O₂

纯度: 98.93%

InChi Key: InChIKey=XWEOIAMCTHLJJB-UHFFFAOYSA-M

Smiles: [I-].C(#CC[N+])(C)(C)COC1=NOCC1

外观: 固体粉末

作用通路: Neuroscience

作用靶点: AChR

溶解性: DMSO: 56.7mg/mL (174.9mM; ultrasonic and warming and heat to 60°C)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Iperoxo is an agonist of Muscarinic AChR and can be used to probe activation-related conformational transitions of muscarinic receptors.