

---

## Lansoprazole sulfone

产品编号: D53947

CAS: 131926-99-3

分子式: C<sub>16</sub>H<sub>14</sub>F<sub>3</sub>N<sub>3</sub>O<sub>3</sub>S

纯度: 100%

InChi Key: InChIKey=TVMJMCGRSSSSDJ-UHFFFAOYSA-N

Smiles: Cc1c(OCC(F)(F)F)ccnc1CS(=O)(=O)c1nc2ccccc2[nH]1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 3.85 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Lansoprazole sulfone is a metabolite of the gastric pump inhibitor Lansoprazole.