

---

## GPR35 agonist 2

产品编号: D53950

CAS: 494191-73-0

分子式: C<sub>17</sub>H<sub>11</sub>FN<sub>2</sub>O<sub>3</sub>S

纯度: 97.18%

InChi Key: InChIKey=SYCKPHBALHXMIR-UHFFFAOYSA-N

Smiles: O=C(O)C1=CC=C(C=C1)C=C2SC(=NC2=O)NC=3C=CC=CC3F

外观: 固体粉末

作用通路: Endocrinology/Hormones;GPCR/G Protein

作用靶点: GPR

溶解性: DMSO : 22.5 mg/mL (65.7 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GPR35 agonist 2 is a potent GPR35 agonist. The EC<sub>50</sub> values of GPR35 agonist 2 in the β-arrestin and Ca<sup>2+</sup> release assays were 26 and 3.2 nM, respectively.