
DLPC

产品编号: D53966

CAS: 18194-25-7

分子式: C₃₂H₆₄NO₈P

纯度: 98.74%

InChi Key: InChIKey=IJFVSSZAOYLHEE-SSEXGKCCSA-N

Smiles: CCCCCCCCCCCC(OC[C@@H](OC(CCCCCCCCCC)=O)COP(OCC[N+](C)(C)C)([O-])=O)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO : 7.14mg/ml (11.48mM) H₂O : <1 mg/ml; Ethanol : 99 mg/mL (159.20 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: DLPC is an LRH-1 agonist ligand. DLPC is a phospholipid for biological study[1][2].