
Bevurogant

产品编号: D53978

CAS: 1817773-66-2

分子式: C₂₆H₂₈N₈O₃S

纯度: 98.78%

InChi Key: InChIKey=HVWHIBHBDCYLDI-HNNXBMFYSA-N

Smiles: [C@@H](C)(N1C=2C(N=C(NCC3=CC=C(S(C)(=O)=O)C=N3)C1=O)=CN=C(N2)C=4C(=NC=NC4C)C5CC5)C6CC6

外观: 固体粉末

作用通路: Metabolism

作用靶点: ROR

溶解性: DMSO : 245 mg/mL (459.99 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bevurogant is an antagonist of ROR γ t receptor and can be used in studies about the treatment of chronic inflammatory diseases.