
(+)-PD 128907 hydrochloride

产品编号: D54001

CAS: 300576-59-4

分子式: C₁₄H₂₀ClNO₃

纯度: 98.71%

InChi Key: InChIKey=DCFXOTRONMKUJB-QMDUSEKHSA-N

Smiles: CCCN1CCO[C@@H]2c3cc(O)ccc3OC[C@H]21.Cl

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Dopamine Receptor

溶解性: DMSO:150 mg/mL (524.90 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: (+)-PD 128907 hydrochloride is an agonist of D3 dopamine receptor .