

---

## L-168049

产品编号: D54024

CAS: 191034-25-0

分子式: C<sub>24</sub>H<sub>20</sub>BrClN<sub>2</sub>O

纯度: 100%

InChi Key: InChIKey=HHBOWXZOLYQFNY-UHFFFAOYSA-N

Smiles: ClC=1C=CC(=CC1)C2=CC(C=3C=C(Br)C=CC3OCCC)=C(N2)C=4C=CN=CC4

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: Glucagon Receptor

溶解性: DMSO: insoluble; Ethanol: 3.4mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: L-168049 is a selective and non-competitive antagonist of human glucagon receptor with IC<sub>50</sub>s of 3.7 nM, 63 nM, and 60 nM for human, murine, and canine, respectively.