
CB-5339

产品编号: D54025

CAS: 1863952-15-1

分子式: C₂₄H₂₄N₆O

纯度: 98.5%

InChi Key: InChIKey=XDHFSLWWYBVSLN-UHFFFAOYSA-N

Smiles: O=C(N)C1=CC=CC2=C1C=C(N2C3=NC(NCC=4C=CC=CC4)=C5C(=NCCC5)N3)C

外观: 固体粉末

作用通路: Ubiquitination

作用靶点: p97

溶解性: DMSO: 37.8mg/mL (91.6mM; ultrasonic and warming and heat to 60°C)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CB-5339 is a potent p97 inhibitor with an IC₅₀ <30 nM (WO2015109285A1, compound FF07).