
CD437

产品编号: D54046

CAS: 125316-60-1

分子式: C₂₇H₂₆O₃

纯度: 98.79%

InChi Key: InChIKey=LDGIHZJOIQSHPB-UHFFFAOYSA-N

Smiles: O=C(C1=CC=C2C=C(C3=CC=C(O)C(C4(C5)CC6CC5CC(C6)C4)=C3)C=CC2=C1)O

外观: 固体粉末

作用通路: Autophagy; Metabolism

作用靶点: Autophagy; Retinoid Receptor

溶解性: DMSO:150 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CD437 is a specific Retinoic Acid Receptor γ (RAR γ) agonist.