
Ibrutinib deacryloylpiperidine

产品编号: D54075

CAS: 330786-24-8

分子式: C₁₇H₁₃N₅O

纯度: 99.8%

InChi Key: InChIKey=XBBITICRHQJPNE-KRWDZBQOSA-N

Smiles: Nc3ncnc4nnc(c2ccc(Oc1ccccc1)cc2)c34

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: BTK

溶解性: DMSO: 10 mg/mL (32.97 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ibrutinib deacryloylpiperidine is a highly selective Bruton' s tyrosine kinase (BTK) irreversible inhibitor with an IC₅₀ of 0.5 nM.