
GSK3987

产品编号: D54090

CAS: 264206-85-1

分子式: C₂₄H₂₀N₂O₃

纯度: 99.48%

InChi Key: InChiKey=HLZMYWLMBBLASX-UHFFFAOYSA-N

Smiles: COc(cc1)ccc1NC(C(N(Cc1ccccc1)C1=O)=O)=C1c1ccccc1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Liver X Receptor

溶解性: DMSO: 50 mg/ml (130.06 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK3987 is an LXR ligand. GSK3987 recruits the steroid receptor coactivator-1 to human LXRA and LXRbeta with EC50s of 40 nM.