
Exarafenib

产品编号: D54294

CAS: 2639957-39-2

分子式: C₂₆H₃₄F₃N₅O₃

纯度: 99.26%

InChi Key: InChiKey=GZMYLSJUNSCMTD-MOPGFXCFSA-N

Smiles: CC=1C(C=2C=C(N=C(N[C@@H](CO)C)C2)N3CCOCC3)=CC(NC(=O)N4C[C@H](CC(F)(F)F)CC4)=CC1

外观: 固体粉末

作用通路: MAPK

作用靶点: p38 MAPK;Raf

溶解性: DMSO: 50 mg/mL (95.86 mM,Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Exarafenib is a potent and orally available pan-RAF inhibitor.Exarafenib has antitumor activity and acts by inhibiting downstream MAPK pathway signaling.Exarafenib is used in cancer research.