
GSK199

产品编号: D54310

CAS: 1549811-53-1

分子式: C₂₄H₂₉ClN₆O₂

纯度: 98.09%

InChi Key: InChIKey=KRGMIOKDGHBYQE-UNTBKODSA-N

Smiles: C(C)N1C(=CC=2C1=NC=CC2)C3=NC=4C(N3C)=C(OC)C=C(C(=O)N5C[C@H](N)CCC5)C4.Cl

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: PAD

溶解性: DMSO: 11 mg/mL (23.45 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK199 is a selective PAD4 inhibitor(IC₅₀ of 200 nM in the absence of calcium).