
BpV(HOpic)

产品编号: D54369

CAS: 722494-26-0

分子式: C₆H₄K₂N₂O₈V

纯度: 99.86%

InChi Key: InChIKey=UBYMWHIAIGQCOR-UHFFFAOYSA-M

Smiles: [K]O[V]12(OO1)(OO2)N3=CC(O)=C/C=C3/C(=O)O[K]

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PTEN

溶解性: DMSO: Insoluble; Water: 10 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: bpV (HOpic) is a potent inhibitor of PTEN (IC₅₀ of 14 nM).