
TH-263

产品编号: D54481

CAS: 313520-94-4

分子式: C₂₁H₂₀N₂O₃S

纯度: 99.33%

InChi Key: InChIKey=QDGVJMITKNOVTP-UHFFFAOYSA-N

Smiles: O=C(NCc1ccccc1)c1ccc(cc1)S(=O)(=O)NCc1ccccc1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 25 mg/mL (65.71 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TH-263 is inactive toward both LIMK1 and LIMK2 and thus can be used as negative control for TH-257, and is a diaryl sulfonamide derivative.