
PKI-402

产品编号: D54488

CAS: 1173204-81-3

分子式: C₂₉H₃₄N₁₀O₃

纯度: 100%

InChi Key: InChIKey=ZAXFYGBKZSQBIV-UHFFFAOYSA-N

Smiles: O=C(NC=1C=CC(=CC1)C=2N=C3C(N=NN3CC)=C(N2)N4CCOCC4)NC5=CC=C(C=C5)C(=O)N6CCN(C)CC6

外观: 固体粉末

作用通路: PI3K/Akt/mTOR

作用靶点: PI3K; mTOR

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PKI-402 is a potent PI3K and mTOR inhibitor. PKI-402 inhibits PI3K α , mTOR, PI3K β , PI3K δ and PI3K γ with IC₅₀s of 2, 3, 7, 14 and 16 nM, respectively.