
DG051

产品编号: D54523

CAS: 929915-58-2

分子式: C₂₁H₂₅Cl₂NO₄

纯度: 98.12%

InChi Key: InChIKey=UCPVOTSNNAVKNE-LMOVXPDSA-N

Smiles: Clc1ccc(Oc2ccc(OC[C@@H]3CCCN3CCCC(O)=O)cc2)cc1.Cl

外观: 固体粉末

作用通路: Metabolism

作用靶点: Aminopeptidase

溶解性: DMSO:317 mg/mL (743.56 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: DG051 is a potent leukotriene A₄ hydrolase (LTA₄H) inhibitor (IC₅₀: 47 nM).