
SR9238

产品编号: D54529

CAS: 1416153-62-2

分子式: C₃₁H₃₃NO₇S₂

纯度: 99.99%

InChi Key: InChiKey=HDZWHJYZJWLTAG-UHFFFAOYSA-N

Smiles: CCOC(c1ccc(CN(Cc(cc2)ccc2-c2cc(S(C)(=O)=O)ccc2)S(c2c(C)cc(C)cc2C)(=O)=O)o1)=O

外观: 固体粉末

作用通路: Metabolism

作用靶点: Liver X Receptor

溶解性: DMSO : 50 mg/mL (83.93 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SR9238 is a synthetic agonist of liver X receptor inverse (IC₅₀s: 214 nM and 43 nM for LXR α and LXR β , respectively).