
LDN-27219

产品编号: D54533

CAS: 312946-37-5

分子式: C₂₀H₁₆N₄O₂S₂

纯度: 99.15%

InChi Key: InChIKey=WLBUIQBNZXIDJ-UHFFFAOYSA-N

Smiles: NNC(=O)CSc1nc2scc(-c3ccccc3)c2c(=O)n1-c1ccccc1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 100 mg/mL (244.8 mM), Need ultrasonic

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LDN-27219 is a potent hTGase inhibitor with IC₅₀ of 0.6 μM.