

---

## C12-200

产品编号: D54547

CAS: 1220890-25-4

分子式: C<sub>70</sub>H<sub>145</sub>N<sub>5</sub>O<sub>5</sub>

纯度: 99.51%

InChi Key: InChIKey=RVHYPUORVDKRTM-UHFFFAOYSA-N

Smiles: OC(CN(CCN1CCN(CC1)CCN(CC(O)CCCCCCCCC)CC(O)CCCCCCCCC)CCN(CC(O)CCCCCCCCC)CC(O)CCCCCCCCC)CCCCCCCCC

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 90 mg/mL (79.2 mM; Need ultrasonic )

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: C12-200 is a benchmark ionizable cationic lipidoid along with helper lipids and can be used in the delivery of mRNA.