
GSK-A1

产品编号: D54549

CAS: 1416334-69-4

分子式: C₂₉H₂₇FN₆O₄S

纯度: 98.97%

InChi Key: InChiKey=AJOGHKUZDLYXKS-UHFFFAOYSA-N

Smiles: COC(ncc(-c(cc1)cc2c1nc(N)n2-c(cc1)ccc1N1CCOCC1)c1)c1S(Nc(ccc1)c1F)(=O)=O

外观: 固体粉末

作用通路: Microbiology/virology; Proteases/Proteasome; PI3K/Akt/mTOR signaling

作用靶点: HCV Protease; PI4K

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK-A1 is a specific inhibitor of PI4KA, also inhibited calcium-driven PRL secretion without affecting calcium signaling and Prl expression. GSK-A1 has anti-hepatitis C virus (HCV) potential.