
PI3Kδ-IN-3

产品编号: D54572

CAS: 1431540-99-6

分子式: C₂₈H₂₄N₆O₃

纯度: 98.43%

InChi Key: InChIKey=QYTXJLQB SYAMGR-UHFFFAOYSA-N

Smiles: N#CC=1C=C(C=NC1OC)C2=CC=C3N=CN=C(C=4C=CC=C(C4)C(=O)N5CCN(C(=O)C)CC5)C3=C2

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PI3Kδ-IN-3 is a PI3Kδ inhibitor with an IC₅₀ value of 9 nM. PI3Kδ-IN-3 has a favorable pharmacokinetic profile and inhibits B cell function.