

---

## PK095

产品编号: D54589

CAS: 380314-37-4

分子式: C<sub>20</sub>H<sub>18</sub>N<sub>4</sub>O<sub>2</sub>S

纯度: 96.84%

InChi Key: InChIKey=ALEUTTNIZDFZSP-UHFFFAOYSA-N

Smiles: O=C1C=CN=C(SCC(=O)NC2=CC=C3C(=C2)C=4C=CC=CC4N3CC)N1

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: ATPase

溶解性: DMSO: 3.79 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PK095 is a proprietary compound in the guanidine - based F1F0-ATPase inhibitor family.