
Tabernanthalog

产品编号: D54592

CAS: 2483829-58-7

分子式: C₁₄H₁₈N₂O

纯度: 98.57%

InChi Key: InChIKey=FNGNYGCPNKZYOG-UHFFFAOYSA-N

Smiles: O(C=1C=CC2=C(C1)NC3=C2CCN(C)CC3)C

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 45 mg/mL (195.39 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tabernanthalog is an analog of ibogaine. Tabernanthalog is a 5-HT_{2A} agonist.