

---

## Mebrofenin

产品编号: D54609

CAS: 78266-06-5

分子式: C<sub>15</sub>H<sub>19</sub>BrN<sub>2</sub>O<sub>5</sub>

纯度: 99.29%

InChi Key: InChIKey=MHPZZZZLAQGTHT-UHFFFAOYSA-N

Smiles: Cc1cc(C)c(NC(=O)CN(CC(O)=O)CC(O)=O)c(C)c1Br

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:250 mg/mL (645.61 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Mebrofenin is available as a ready to use the kit for radio-labeling with Tc-99m. Technetium Tc 99m Mebrofenin is a diagnostic radiopharmaceutical composed of diisopropyl-iminodiacetic acid (DISIDA) attached to a technetium-99m ion.