
PD128907 HCl

产品编号: D54621

CAS: 112960-16-4

分子式: C₁₄H₂₀ClNO₃

纯度: 98%

InChi Key: InChIKey=DCFXOTRONMKUJB-KYSPHBLOSA-N

Smiles: C(CC)N1[C@@]2([C@](C=3C(OC2)=CC=C(O)C3)(OCC1)[H])[H].Cl

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Dopamine Receptor

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PD128907 HCl is a potent and selective dopamine D3 receptor agonist, with EC₅₀ of 0.64 nM, exhibits 53-fold selectivity over dopamine D2 receptor.