
STING (WT/R232) activity regulator 126

产品编号: D54634

CAS: 702662-64-4

分子式: C₂₃H₁₈ClFN₂O₂S

纯度: 99.61%

InChi Key: InChiKey=VJXUTYSJHGFAKG-UHFFFAOYSA-N

Smiles: O=C(c(cc1)cc(N2Cc(c(F)ccc3)c3Cl)c1SCC2=O)NCc1cccc1

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: STING

溶解性: DMSO : 55.00mg/ml

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: STING (WT/R232) activity regulator 126 is a STING modulator with antibacterial and antiviral effects.