

---

## AC-262536

产品编号: D54678

CAS: 870888-46-3

分子式: C<sub>18</sub>H<sub>18</sub>N<sub>2</sub>O

纯度: 99.38%

InChi Key: InChIKey=ATKWLNSCJYLXPF-FICVDOATSA-N

Smiles: N#CC1=C2C=CC=CC2=C(C=C1)N3[C@H]4C[C@@H](C[C@@H]3CC4)O

外观: 固体粉末

作用通路: Endocrinology/Hormones

作用靶点: Androgen Receptor

溶解性: DMSO : 50 mg/mL (179.63 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AC-262536 is a selective androgen receptor (SAR) modulator and exhibits potent agonist activity at the androgen receptor. It has partial agonist activity with respect to testosterone and suppresses luteinizing hormone.