
ML604440

产品编号: D54713

CAS: 1140517-08-3

分子式: C₁₇H₂₄BF₃N₂O₄

纯度: 95.24%

InChi Key: InChIKey=IRZZVDHGNHGNAL-ZDUSSCGKSA-N

Smiles: OB([C@@H])(NC(C(NC(C1=CC=CC=C1C(F)(F)F)=O)(C)C)=O)CC(C)C)O

外观: 固体粉末

作用通路: Proteases/Proteasome; Ubiquitination

作用靶点: Proteasome

溶解性: DMSO:100 mg/mL (257.61 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ML604440 is a cell permeable proteasome β 1i (LMP2) subunit inhibitor.