
FAUC 213

产品编号: D54722

CAS: 337972-47-1

分子式: C₁₈H₁₉ClN₄

纯度: 99.78%

InChi Key: InChIKey=DTRXURJDKOYCCD-UHFFFAOYSA-N

Smiles: C1CN(CCN1CC2=NN3C=CC=CC3=C2)C4=CC=C(C=C4)Cl

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Dopamine Receptor

溶解性: DMSO: 3.26 mg/mL (9.99 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: FAUC-213 is a selective full antagonist of the dopamine D4 receptor.