

---

## CU-115

产品编号: D54727

CAS: 2471982-20-2

分子式: C<sub>21</sub>H<sub>11</sub>F<sub>7</sub>INO<sub>2</sub>

纯度: 100%

InChi Key: InChIKey=QYSLCRYVUJORPX-UHFFFAOYSA-N

Smiles: O=C(NC1=CC=C(OC=2C=C(C=C(C2)C(F)(F)F)C(F)(F)F)C=C1)C=3C(F)=CC=CC3I

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: TLR

溶解性: DMSO : 90 mg/mL (158.1mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CU-115 is a selective antagonist of TLR8 with IC<sub>50</sub>s of 1.04 μM and >50 μM for TLR8 and TLR7, respectively.