

**bs-10499R**

**[ Primary Antibody ]**

## Listeriolysin Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Listeriolysin

**Immunogen:** KLH conjugated synthetic peptide derived from Listeria monocytogenes Listeriolysin: 251-350/529.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** **IHC-P** (1:100-500)

**IHC-F** (1:100-500)

**IF** (1:100-500)

**ICC/IF** (1:100-500)

**ELISA** (1:5000-10000)

**Reactivity:** (predicted: Listeria monocytogenes)

**Predicted MW.:** 55 kDa

**Subcellular Location:** Secreted ,Cell membrane