
Qc1

产品编号: D54783

CAS: 403718-45-6

分子式: C₂₃H₁₆F₃N₃O₂S

纯度: 100%

InChi Key: InChIKey=IFNVTSPDMUUAIFY-UHFFFAOYSA-N

Smiles: O=C(NCC=1C=CC=CC1)C2=CC=C3C(=O)N(C4=CC=CC(=C4)C(F)(F)F)C(=S)NC3=C2

外观: 固体粉末

作用通路: Metabolism

作用靶点: Dehydrogenase

溶解性: DMSO: 20 mg/ml (43.9 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Qc1 is an inhibitor of threonine dehydrogenase (TDH) inhibitor and can be used in studies about metabolic diseases.