
DL-AP5

产品编号: D54798

CAS: 76326-31-3

分子式: C₅H₁₂N₂O₅P

纯度: 98%

InChi Key: InChIKey=VOROEQBFPPPIACJ-UHFFFAOYSA-N

Smiles: NC(CCCP(O)(O)=O)C(O)=O

外观: 固体粉末

作用通路: Neuroscience

作用靶点: NMDAR

溶解性: DMSO:Insoluble;Water:10 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: DL-AP5 is the racemic version of the selective N-methyl-D-aspartate (NMDA) receptor antagonist, with anticonvulsant.