
Larotinib

产品编号: D54829

CAS: 1438072-11-7

分子式: C₂₄H₂₆ClFN₄O₄

纯度: 98.45%

InChi Key: InChIKey=UFUKZMIXTOMMLU-ZRZAMGCNSA-N

Smiles: N(C=1C2=C(C=C(OC)C(OCCCN3C[C@@]4([C@](C3)(OCCO4)[H])[H])=C2)N=CN1)C5=CC(Cl)=C(F)C=C5

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors;Angiogenesis;JAK/STAT signaling;Tyrosine Kinase/Adaptors;Angiogenesis;Tyrosine Kinase/Adaptors;Immunology/Inflammation;NF-Kb

作用靶点: Tyrosine Kinases;EGFR;BTK;IRAK

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Larotinib is an orally active, potent, and broad-spectrum tyrosine kinase inhibitor (TKI) with an IC₅₀ of 0.6 nM for EGFR.