
KCNQ2/3 activator-1

产品编号: D54846

CAS: 1009344-33-5

分子式: C₂₃H₂₉N₂O

纯度: 98.99%

InChi Key: InChIKey=FJNPZKZPWWSON-UHFFFAOYSA-N

Smiles: CC(C)(C)CC(Nc(c(C)c1)c(C)cc1N(CC1)Cc(cc2)c1cc2F)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 45 mg/mL (122.12 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: KCNQ2/3 activator-1 is a Kv7.2/Kv7.3 (KCNQ2/3) potassium channel activator, details from Compound A of patent WO2021113757A1.