
PU139

产品编号: D54852

CAS: 158093-65-3

分子式: C₁₂H₇FN₂OS

纯度: 97.61%

InChi Key: InChIKey=QMCIVCACYJRAAY-UHFFFAOYSA-N

Smiles: O=C1C2=CC=CN=C2SN1C3=CC=C(F)C=C3

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Acetyltransferase

溶解性: DMSO : 10mg/mL (40.61 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PU139 is a novel inhibitor of histone acetyltransferase (HAT).