
SAR405 R enantiomer

产品编号: D54858

CAS: 1946010-79-2

分子式: C₁₉H₂₁ClF₃N₅O₂

纯度: 98.04%

InChi Key: InChIKey=SPDQRCUBFSRAFI-IUODEOHRSA-N

Smiles: O=C1C=C(N2[C@H](C)COCC2)N=C3N1CC[C@H](C(F)(F)F)N3CC4=CC(Cl)=CN=C4

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:100 mg/mL (225.30 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SAR405 R enantiomer is the less active enantiomer of SAR405. SAR405 is an inhibitor of PIK3C3/Vps34.