
TC LPA5 4

产品编号: D54888

CAS: 1393814-38-4

分子式: C₂₃H₂₃ClN₂O₃

纯度: 99.05%

InChi Key: InChIKey=BNALUYKEGYUHQ-C-UHFFFAOYSA-N

Smiles: O=C(O)C1=NN(C=CC=C(OC)C2C(=C1)C=CC(=C(Cl)C3)C4CCCCC4

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: LPA Receptor

溶解性: DMSO : 56.2 mg/mL (136.9 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TC LPA5 4 is a LPA5 receptor antagonist.