
Bevenopran

产品编号: D54913

CAS: 676500-67-7

分子式: C₂₀H₂₆N₄O₄

纯度: 98%

InChi Key: InChIKey=ZGCYVRNZWGUXNQ-UHFFFAOYSA-N

Smiles: O=C(N)C=1N=CC(=NC1)OC2=CC=C(C=C2OC)CNCCC3CCOCC3

外观: 固体粉末

作用通路: Endocrinology/Hormones;GPCR/G Protein;Neuroscience

作用靶点: Opioid Receptor

溶解性: DMSO : 120 mg/mL (310.53 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bevenopran is a peripheral antagonist of μ -opioid receptor and can be used in studies about the treatment of opioid-induced bowel dysfunction.