
MS0124

产品编号: D54932

CAS: 1197196-63-6

分子式: C₂₀H₂₉N₅O₃

纯度: 99.56%

InChi Key: InChIKey=MEVMMQQJOXBSAX-UHFFFAOYSA-N

Smiles: N=1C(=NC(NC2CCN(C)CC2)=C3C=C(OC)C(OC)=CC13)N4CCOCC4

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Methyltransferase

溶解性: DMSO: 3.87 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MS0124 is a highly selective G9A-like protein (GLP) inhibitor with an ic₅₀ of 13±4 nM against GLP and 440±63 nM against G9a.