
BMS817378

产品编号: D54969

CAS: 1174161-69-3

分子式: C₂₄H₁₈ClF₂N₄O₇P

纯度: 100%

InChi Key: InChIKey=KXDZWUPUSDCGDD-UHFFFAOYSA-N

Smiles: Nc1nccc(Oc2ccc(NC(=O)c3cn(COP(O)(O)=O)cc(-c4ccc(F)cc4)c3=O)cc2F)c1Cl

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR

溶解性: DMSO: 5.79 mg/ml (10 mM, Need ultrasonic); Water: Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BMS817378 is a potent and selective inhibitor of MET (IC₅₀: 1.7 nM).