
TPPU

产品编号: D55050

CAS: 1222780-33-7

分子式: C₁₆H₂₀F₃N₃O₃

纯度: 98%

InChi Key: InChIKey=AAJMQTLFRTZCJK-UHFFFAOYSA-N

Smiles: CCC(=O)N1CCC(CC1)NC(=O)Nc1ccc(OC(F)(F)F)cc1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Epoxide Hydrolase

溶解性: DMSO:28 mg/mL (77.92 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TPPU is a potent inhibitor of both human and mouse sEH (IC₅₀ of 3.7 and 2.8 nM, respectively)