

---

## 1-phenyl-3,4-dihydroisochinoline

产品编号: D55054

CAS: 52250-50-7

分子式: C<sub>15</sub>H<sub>13</sub>N

纯度: 98.01%

InChi Key: InChIKey=CTOQBSUYGFNMJX-UHFFFAOYSA-N

Smiles: C1c(cccc2)c2C(c2ccccc2)=NC1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 65 mg/ml (313.6 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 1-phenyl-3,4-dihydroisochinoline is a structural analog of the neurotransmitter dopamine that interacts with dopamine receptors in the central nervous system. It has been studied as a potential therapeutic agent for the treatment of neurological disorders such as Parkinson's disease, Alzheimer's disease and Huntington's chorea.