
Dilmapimod

产品编号: D55077

CAS: 444606-18-2

分子式: C₂₃H₁₉F₃N₄O₃

纯度: 98.61%

InChi Key: InChIKey=ORVNHOYNEHYKJG-UHFFFAOYSA-N

Smiles: Cc(cc(cc1)F)c1-c1c(C=CC(N2c(c(F)ccc3)c3F)=O)c2nc(NC(CO)CO)n1

外观: 固体粉末

作用通路: Autophagy; MAPK

作用靶点: Autophagy; p38 MAPK

溶解性: DMSO:125 mg/mL (273.87 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Dilmapimod is a potent inhibitor of p38 MAPK ,it potentially suppresses inflammation in chronic obstructive pulmonary disease.