
Locostatin

产品编号: D55093

CAS: 133812-16-5

分子式: C₁₄H₁₅NO₃

纯度: 99.84%

InChi Key: InChIKey=UTZAFVPPWUIPBH-UHFFFAOYSA-N

Smiles: CC=CC(=O)N1C(COC1=O)CC2=CC=CC=C2

外观: 固体粉末

作用通路: MAPK

作用靶点: Raf

溶解性: ethanol:10 mg/ml;DMSO:5 mg/ml;DMF:5 mg/ml

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Locostatin is a potent and cell permeable inhibitor of Raf kinase inhibitor protein (RKIP)/Raf1 kinase interaction and an inhibitor of cell migration.