
BCI-137

产品编号: D55123

CAS: 112170-24-8

分子式: C₁₁H₁₁N₃O₆S

纯度: 95%

InChi Key: InChIKey=VMSRTLYZOWBBFW-UHFFFAOYSA-N

Smiles: CC(NS(=O)(=O)c1ccc2[nH]c(=O)c(=O)[nH]c2c1)C(O)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 65 mg/ml (207.48 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BCI-137 is an Argonaute 2 ligand. It acts by targeting the MicroRNA Binding Domain.