

---

## AR420626

产品编号: D55130

CAS: 1798310-55-0

分子式: C<sub>21</sub>H<sub>18</sub>Cl<sub>2</sub>N<sub>2</sub>O<sub>3</sub>

纯度: 98%

InChi Key: InChIKey=GGTYQECCGLBHGS-UHFFFAOYSA-N

Smiles: CC(NC(CCC1)=C(C2c3ccco3)C1=O)=C2C(Nc(cc(cc1)Cl)c1Cl)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 10mg/mL (Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AR420626 is a selective fatty acid 3 receptor (FA3R) agonist and provides protective effects against salsolionl (SALS), which was totally blocked by beta-hydroxy butyrate (BHB), a selective FA3R antagonist.